

Device	Company	Accuracy MARD (Low is better)	Needs Fingertick Calibration	Alarms for Lows	Connects with a Pump	Medicare Approval	Days of Use Per Sensor	Connect to Phone	Comment
Dexcom G6	Dexcom	9.3%	No	Yes	Displays numbers on Tandem and OmniPod pumps. Can prevent lows on latest Tandem pump.	No	10 days	Yes	Awaiting Medicare approval.
Dexcom G5	Dexcom	9.0%	2 times a day	Yes	Displays numbers on Tandem and OmniPod pumps.	Yes	7 days	Yes	
FreeStyle Libre	Abbot	9.4%	No	No	No	Yes	14 days	No	Cash discounts available. A site is: https://www.singlecare.com/
Eversense	Sensionics	8.5%	2 times a day	Yes	No	No	3 months	Yes	Implanted under the skin needs to be replaced every 3 months.
Guardian Connect	Medtronic	9.1%	2-3 times a day	Yes	Yes. With Medtronic 670G provides protecting against lows and adjusts insulin to maintain glucose when not eating. Can display glucose on other pumps.	No	6 days	Yes	Works best with Bayer Contour test strips for calibration. Has Apple phone app to analyze data and give advice using IBM Watson.

If you are considering a CGM, please check for new information. Every effort was made to be accurate as of 10/08/18.