



Dopamine Dose Chart

DOPAMINE DRIP RATES

Dopamine concentration = 1600 mcg/ml solution = 400 mg in 250 ml D5W or NS

Drops per minute based on microdrip tubing (60 gtt/ml)

Patient Weight (kg)	5 mcg / kg / min	10 mcg / kg / min	15 mcg / kg / min	20 mcg / kg / min
40	8	15	23	30
45	8	17	25	34
50	9	19	28	38
55	10	21	31	41
60	11	23	34	45
65	12	24	37	49
70	13	26	39	53
75	14	28	42	56
80	15	30	45	60
85	16	32	48	64
90	17	34	51	68
95	18	36	53	71
100	19	38	56	75
105	20	39	59	79
110	21	41	62	83